

Cordless Fall Monitor



Model #
TL-2100GR1 (1 yr)
\$119.95

The NEW CordLess® Fall Monitor Works with two sensor pads simultaneously! The cordless pad technology communicates with the monitor using wireless technology. Use one CordLess® monitor with two sensor pads for a safer and more complete resident monitoring environment.

Bed and Chair Weight Sensitive Pads

Cordless Chair Pressure Pad



GCT-45 10"x15" (45 uw) - \$61.95
GCT-90 10"x15" (90 uw) - \$69.95
GCT-WI 10"x15" (1 uw) - \$89.95

Cordless Bed Sensor Pad



GBT-45 10"x30" (45 uw) - \$73.95
GBT-90 10"x30" (90 uw) - \$79.95
GBT-RI 10"x30" (1 uw) - \$89.95

Cordless Bed Sensor Pad



GBT-WI 20"x30" (1 uw) \$93.95

When connected to a fall monitor, the pressure pad will trigger the alarm when weight is removed from the pad. These devices are used in long-term care facilities, hospitals, and other situations to help caregivers to monitor residents that may be at risk of falling or wandering. "Chair pads" are pressure pads that are sized for use in chairs and wheelchairs. When a resident gets out of their chair or wheelchair, the chair pad tells the monitor to alarm, thereby notifying a caregiver that the resident is standing up. "Bed pads" are pressure pads that are sized for use on beds. When a resident gets out of bed, the bed pad tells the monitor to alarm, thereby notifying a caregiver that the resident is getting up.

Early Warning Seatbelt

Buckle Seatbelt: TL-2109



Wheelchair seatbelt features a one finger, easy-release buckle connected to a fall monitor. When the belt is unbuckled, the monitor will sound an alert. For use with the TL-2100 series and TL-2016 wireless monitors.

Adjustable belt: 20" to 54"
\$32.95 – 1yr warranty

Hook and Loop Seatbelt: TL-2109V



Wheelchair seatbelt features a Hook and Loop, quick release belt connected to a fall monitor. When the belt is unfastened, the monitor will sound an alert. For use with the TL-2100 series and TL-2016 wireless monitors.

Seat Belts are Non-restraining.
\$35.95 – 1yr warranty

Pressure-Sensing Wheelchair

The Wheelchair Straight Cushion reduces forward slipping and helps the resident to have more comfortable posture. For added security, this is a pressure-sensing cushion. The cushion is to be used with pressure pad monitors (sold separately) to notify the caregiver when the resident attempts to vacate his/her chair.



CC-01-\$59.95
Straight cushion-30 day warranty

CC-02-\$107.95
Pressure-Sensing w/alarm function-30 day warranty

Fall Mattress

Cordless Pressure-Sensitive Fall Mattress w/Alarm Function

This fall mattress is particularly useful to protect the residents from fall-related injuries in facilities that prohibit the use of side rails on beds. This firm yet soft, 2" thick mattress spreads the weight from impact, helping to protect knees, hips and elbows. Designed to reduce the chance of serious injury. This mattress can be used with any FallGuard® monitor to notify the caregiver if a resident falls. Folding design allows for easy storage. 30 day warranty.



Cordless Pressure-Sensing w/Alarm Function – GWC-03 – (2 sections 36"x68"x2") - \$239.95

Pressure-Sensing w/Alarm Function – WC-03 –(2 sections 36"x68"x2) - \$163.95

Fall Mattress – WC-01 –(3 sections) -\$83.95

Weight-Sensitive Floor Mat for Standard or Wireless Fall Prevention/Mobility Alarm Monitor

A weight-sensitive floor mat connects to an alarming monitor and sets off the alarm when weight is placed on the floor mat. The floor mat may be placed at the side of a bed or in a doorway to alert the caregiver when a resident arises or attempts to depart. Our floor mats feature phone-style plugs. Simply plug the floor mat into a compatible monitor and it is ready for use.



\$119.95-1 year warranty
Standard or wireless (FM-03 black)24"x36"

\$167.95- 1 year warranty
Wireless (FMT-07C Grey) 24"x48"

\$119.95-1 year warranty
Standard or wireless (FM-05 Gray)24"x36"

\$119.95 – 1 year warranty
Extra- long safe TRelease Cord (FM-09Gray) 24"x64"

\$119.95- 1 year warranty
Standard or wireless safe TRelease (FM-07Gray)24"x48"

\$89.95 – 6 month warranty
Standard or wireless (FM-06grey/vinyl)20"x30"

Fall Monitors



Cordless® Fall Monitor

Safety auto-reset™, status light, low battery and pad lost alerts, adjustable tone, volume, and delay.
For use with Cordless Bed and Chair Sensor Pads and Floor Mats

1 year warranty
TL-2100G - \$119.95



Recordable Voice Fall Monitor

Safety auto-reset™, status light, low battery, and pad lost alerts, adjustable tone, volume, and delay.
Use with Corded Bed Sensor Pads, Chair Sensor Pads, Floor Mats and Seat Belts

1 year warranty
TL-2100V - \$83.95



Safety Auto-Reset™ Fall Monitor

Status light, low battery, and lost pad alerts, adjustable tone, volume, and delay. For Corded Bed Sensor Pads, Chair Sensor Pads, Floor Mats and Seat Belt

1 year warranty
TL-2100S - \$71.95



Economy Fall Monitor

On/off switch, pad lost and low battery indicator, adjustable tone and volume. For Corded Bed Sensor Pads, Chair Sensor Pads, Floor Mats and Seat Belts

1 year warranty
TL-2100E - \$59.95

Change Pad Monitor

No more Dating Pads! Helps Avoid F-Tags!! TL-2100CP - \$119.95

The NEW Change Pad Indicator™ Monitor alerts the caregiver when a Timed™ pressure pad has expired. The federal government requires care facilities to only use pressure pads that are within their warranty period. Generally, caregivers write the date of purchase on the pad with a marker. They then have to continually check the dates on all the pads in their facility and keep up with which pads need to be changed and when. The Change Pad Indicator™ Monitor, when used with a Timed™ pressure pad, will keep track of the age of the pad and will indicate (with a light and a sound) when the pad has expired. This is a tremendous help to caregivers because government inspectors will write up facilities (giving them “F-Tags”) if pressure pads in use are past their warranties. The Timed™ pressure pads for bed or chair function like normal pressure pads, but the monitor will alert the caregiver when it is time to change the pads. Starting 2 weeks before the pad expires, the “Change Pad” light will blink and the alarm will beep to indicate that the pad’s warranty is almost expired. During this 2 week period, the caregiver can reset (silence) the alarm and then the alarm and light will activate again 24 hours later. 2 weeks later, when the pad has actually expired, the light will stay on and the monitor will alarm, but the caregiver will be unable to silence it. Timed™ pads will work with other monitors, but only the TL-2100CP monitor features a “change pad” indicator.



Timed™ Pressured Sensor pads

“Timed” Bed Sensor Pad SafeTRelease Cord



TB-45R 10”x30” (45 uw) - \$33.95
TB-90R 10”x30” (90 uw) - \$39.95
TB-RI 10”x30” (1 uw) - \$55.95

“Timed” Bed Sensor Pad SafeTRelease Cord



TB-45W 20”x30” (45 uw) - \$33.95
TB-90W 20”x30” (90 uw) - \$43.95
TB-WI 20”x30” (1 uw) - \$59.95

“Timed” Chair Sensor Pad Coiled Cord



GCT-45 10”x15” (45 uw) - \$21.95
GCT-90 10”x15” (90 uw) - \$33.95
GCT-WI 10”x15” (1 uw) - \$47.95

Tidy Toilet Pressure Pads

For use with the TL-2100 Series Monitors
PPT-1 (1 pad) - \$17 PPT-5 (5 pack) - \$87.95

This Pressure Pad is used with TL-2100 Series Fall Monitors. The Fall Alarm System alerts the Caregiver when the resident has gotten off of the toilet, allowing resident more privacy using the restroom. Toilet Seat Pads are placed underneath the toilet seat and when pressure is released off of the toilet seat, the fall alarm sounds. The Toilet Pressure Pads each have a 30 Day Warranty.



Economy Central Monitoring Unit (CMU)

(1-30 device capability)
433MHz – 1 year warranty
433-CMU - \$119.95




The Economy Central Monitor System is programmable for many residents at one central location and can be activated by a variety of peripheral wireless devices so that the caregiver can meet each resident's particular needs. With the use of the 433-NCA adapter, this system can also be used with ordinary pressure pads, weight-sensitive floor mats, and early warning seat belts. The 433-CMU Economy Central Monitor System alarms audibly and the display indicate which resident originated the alarm. The caregiver can be notified by a Caregiver Pager as long as they are located within 200 feet (or more with a bumper) of the central monitor.

Wireless Reset Button



433MHz – 90 day warranty
433-RB - \$33.95

Nurse Call Button



433MHz – 90 day warranty
433-NC - \$21.95

Local Pocket Pager w/ LCD Display



433MHz – 6 month warranty
433-PGD - \$35.95



Window/Door Exit Transmitter

433MHz – 1 year warranty
433-EXT - \$37.95



Signal Bumper between Economy CMU & Pager

Increases range up to 100 ft
1 year warranty
433-BMPG - \$139.95



Motion Sensor
433MHz – 90 day warranty
433-MS - \$35.95



Motion Sensor and Paging System



This simple, streamlined system comes complete with one motion sensor and one Caregiver Pager. The Caregiver Pager has an audible alert and is small enough to be carried in a pocket or clipped to a belt. This Motion Sensor Paging System alerts the caregiver wearing the pager of a resident's activities when the resident interrupts the motion sensor's infrared field, as long as they are within 300 feet of each other. The Motion Sensor may be mounted next to a bed to indicate when a resident arises, or by a doorway to indicate when a resident is approaching the door way.

TL-5102MP - 90 day warranty - \$59.95

Two Call Button Paging System



This simple, streamlined system comes complete with two nurse call buttons and one caregiver pager. The caregiver pager has an audible alarm and is small enough to be carried in a pocket or clipped to a belt. This Two Pendant Paging System alerts the caregiver of a resident's needs when the resident presses one of the nurse call buttons (within 300 feet). The buttons may be mounted on the headboard of a bed or worn as a pendant around the resident's neck. The caregiver paging system can also be configured for use with chair and bed pressure pads, early warning seatbelts, and pressure-sensitive floor mats. The pager requires 2 AA batteries. The nurse call buttons each require 1 A23-12V. Protective Pager Case Included!

TL-5102TP - 90 day warranty - \$59.95

Bed Safety Wedge



WC-04 - \$47.95
30 day warranty

Wheelchair Straight Cushion



CC-01 - \$59.95
30 day warranty

Fall Mattress – 3 sections



WC-01 - 36" x 68"x2" - \$83.95
30 day warranty

Accessories

Caregiver Key™



Replacement/Additional Caregiver Key™ for the FallGuard® series alarm monitors. When the fall prevention or mobility monitor is set for its use, only the Caregiver Key™ may silence the alarm.

Use with: TL-2100S, TL-2100CP,
TL-2100V
CK-01 - \$2.49

AC Adapter
for Bed and Chair Alarms



Used with any of the bed alarms, chair alarms, door alarms, or wireless bed alarms, the AC-02 eliminates the need for replacing batteries and is UL approved.

Use with: TL-2100E, TL-2100CP,
TL-2100S, TL-2100V
AC-02 - \$9.95

Steel Mounting Bracket Adjustable
for bed or chair



Easily hooks onto Protective Boot for Fall Monitors (included with monitor) and mounts monitor to headboard and rails. Makes the alarm unit easy to move and easy to mount

Use with: TL-2100 series Series
Monitor for Headboards and Rails
MB-03 - \$21.95
MB-03B

Mounting Bracket
for wheelchairs



Use with: TL-2100 series
MB-07B - \$9.95



To order call 510-999-6969

Please allow 5 – 10
business days for delivery

Note : *dw – days warranty
*yw – year warranty

Sales Department
1-877-877-0197
Customer Service
1-800-650-3637
Technical Support
1-800-650-3637

Mon-Fri: 7:30am – 3:00pm PST

Mounting Bracket
for head boards and bed rails



Use with: TL-2100 series
MB-07C - \$9.95